



WE LISTEN.
WE THINK DIFFERENTLY.
WE FIND SOLUTIONS.

Medical devices that speak the language of the users.



ETTFIX ONE or ETTFIX SILIKON



- One-piece strap for fixing endotracheal tubes
- Fixation: placement of the strap on the tube or at the back of the head – Velcro strip must always stick to the tube!
- ETTFix Silicone:
 - Further development of ETTFix One with a continuous silicone layer applied to the elastic band to optimise adhesion



- To prevent unintentional extubation or tube slipping during endotracheal treatments
- ETTFix silicone: particularly suitable for patients with high quantities of saliva due to its strong adhesion

HOW TO SELL

Icebreaker: Questions to ask your customers

- How do you currently fix endotracheal tubes?
- Have there been any skin problems or dislocations?

Target group: This is where ETTFix or ETTFix silicone is needed

Intensive Care Unit | Anaesthesia | A&E | IMC

Benefits: for clinics and clinic staff

- Greater patient safety
- Fewer reintubations
- Time saving

ITEM INFORMATION

Description	note	item number	packaging unit (PU)
ETTFix One	one-piece, white	HSOETT2000	1 PU = 25 pieces
ETTFix Silikon	one-piece, blue	HSOETT2010	1 PU = 25 pieces

**WE ARE THE HOME OF HEALTHCARE.
WE ARE THE HEINTEL GROUP.**

