




WE LISTEN.  
WE THINK DIFFERENTLY.  
WE FIND SOLUTIONS.

Medical devices that speak the language of the users.



# HSO SafeHeat SET



- SafeHeat Pad: generates constant heat by means of a chemical reaction with air
- SafeHeat glove: optimally distributes heat to the limbs
- Constant heat on the limbs for up to 12 hours – safe and hygienic
- Single-use product



- Bedridden people with a pronounced feeling of cold in the limbs



- Glove is made of cotton fleece
- Pad and its filling is made up of iron, water, diatomite, activated carbon, sodium chloride

## HOW TO SELL

### Icebreaker: Questions to ask your customers

- What do you do if patients have cold hands or feet?
- What solutions are you currently using for centralised patients?
- How do you ensure that the warming of the limbs is controlled and safe?

### Target group: This is where the SafeHeat is needed

Intensive Care Unit | Anaesthesia | Neurology | Oncology

### Benefits: for clinics and clinic staff

- Safe warming of the limbs
- Replacement of any unsafe temporary solutions and the associated time savings

## ITEM INFORMATION

Description	item number	packaging unit (PU)	note
HSO SafeHeat Set	HSOSAFEHESETN	1 PU = 1 packet	sufficient for 20 applications
SafeHeat Glove	HSOWHANDO2	1 PU = 1 packet	sufficient for 20 applications
SafeHeat Pad	HSOWPADO1N	1 PU = 1 packet	sufficient for 20 applications

**WE ARE THE HOME OF HEALTHCARE.  
WE ARE THE HEINTEL GROUP.**

