




WE LISTEN.
WE THINK DIFFERENTLY.
WE FIND SOLUTIONS.

Medical devices that speak the language of the users.



SpO₂ SENSORS



SpO₂ BUTTERFLY SENSOR



SpO₂ SENSOR L-SHAPE



SpO₂ SENSOR FINGERSENSOR



SpO₂ SENSOR SOFT CAP „FROG SENSOR“

HOW TO SELL

Icebreaker: Questions to ask your customers

- Have there been any cable breaks or contact problems?
- How reliable are the measured values in everyday clinical practice?
- Have you had any discussions about the cost per patient or per day?
- Is the handling simple and efficient for the nursing staff?

Target group: This is where the SpO₂ sensors are needed

Intensive Care Unit | IMC | Anaesthesia | Surgery | A&E | Cardiology | Geriatrics

Benefits: for clinics and clinic staff

- Reliable measured values even when moving or in difficult patient conditions
- Proven technical compatibility with existing systems
- Stable measurements

ITEM INFORMATION

| Description | note | item number | packaging unit (PU) |
|-----------------------------------|--|-------------------|---------------------|
| SpO ₂ Butterfly Sensor | suitable for Nellcor-OxiMax™, HP, Mindray and other Nellcor-OxiMax™-compatible modules | HSOA1418-SA03MD | 1 PU = 25 pieces |
| SpO ₂ Butterfly Sensor | suitable for Nellcor-OxiMax™, HP, Mindray and other Nellcor-OxiMax™-compatible modules | HSOA1418-SA03* | 1 PU = 25 pieces |
| SpO ₂ Sensor L-Shape | Nellcor-OxiMax™ | HSOA1418-SN07MD* | 1 PU = 25 pieces |
| SpO ₂ Sensor L-Shape | GE | HSOA2520-SN07 | 1 PU = 25 pieces |
| SpO ₂ Frog sensor | HP/Phillips | HSOA0816-SA106PU* | 1 PU = 1 piece |
| SpO ₂ Frog ensor | Nellcor-OxiMax™ | HSOA1418-SA203MU | 1 PU = 1 piece |
| SpO ₂ Fingersensor | mit GE Stecker | HSOA2520-SA03 | 1 PU = 25 pieces |

*These SpO₂ sensors can be ordered on request.

**WE ARE THE HOME OF HEALTHCARE.
WE ARE THE HEINTEL GROUP.**

